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<b>EN ISO 13485</b> "Medical devices - Quality management systems - Requirements for regulatory purposes"
<b>EN ISO 14971</b> "Medical devices - Application of risk management to medical devices"
<b>EN ISO 15223-1</b> "Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements "
<b>EN ISO 10993-10</b> "Biological evaluation of medical devices Part 10: Tests for irritation and skin sensitization"
<b>EN ISO 10993-12</b> "Biological evaluation of medical devices Part 12: Sample preparation and reference materials"
<b>EN ISO 10993-17</b> "Biological evaluation of medical devices Part 17: Establishment of allowable limits for leachable substances"
<b>EN ISO 10993-18</b> "Biological evaluation of medical devices Part 18: Chemical characterization of materials"
<b>EN ISO 10993-23</b> "Biological evaluation of medical devices Part 23: Tests for irritation"
<b>EN ISO 11137-1</b> "Sterilization of health care products — Radiation — Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices"
<b>EN ISO 11137-2</b> "Sterilization of health care products — Radiation — Part 2: Establishing the sterilization dose"
<b>EN ISO 11737-1</b> "Sterilization of health care products — Microbiological methods — Part 1: Determination of a population of microorganisms on products"
<b>EN ISO 11737-2</b> "Sterilization of health care products — Microbiological methods — Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process"
<b>EN ISO 11607-1</b> "Packaging for terminally sterilized medical devices — Part 1: Requirements for materials, sterile barrier systems and packaging systems"
<b>EN ISO 11607-2</b> "Packaging for terminally sterilized medical devices — Part 2: Validation requirements for forming, sealing and assembly processes"

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**Only to be administered by appropriately trained healthcare professionals who are qualified or accredited in accordance with national law.**

**Algeness®** is an absorbable filler indicated for the correction of tissue and / or hypodermic atrophies. **Algeness®** is an injectable gel consisting of agarose.

**Algeness®** is contained in a special "total emptying" double syringe that allows you to mix the product before use and to use all the gel contained in the syringe.

#### COMPOSITION

**Algeness® LD** contains agarose (1%), water (89,1%) and phosphate buffer (9,9%).

**Algeness® HD** contains agarose (1,5%), water (88,7%) and phosphate buffer (9,8%).

**Algeness® VL** contains agarose (2,5%), non cross-linked hyaluronic acid (0,5%), water (87,3%) and phosphate buffer (9,7%).

**Algeness® DF** contains agarose (3,5%), non cross-linked hyaluronic acid (0,4%), water (86,4%) and phosphate buffer (9,7%).

Agarose with a molecular weight between 123-133 kDa (gel permeation chromatography), chosen for its physical and rheological properties suitable for the application site, is able to form an extrudable gel consisting of a three-dimensional matrix capable of mechanically retaining molecules and organic liquids in dynamic equilibrium with its absorption. Therefore, once the agarose has been injected, thanks to its structure, it is able to give viscoelastic support to the extracellular matrix of the tissue by mechanically filling the spaces involved in subdermal atrophy.

**Algeness®**, injected into the subdermal tissue, provides a viscoelastic supplement to the matrix. Agarose restores the lost tissue volume of both the adipose tissue and the connective stroma, ensuring biocompatibility with the extracellular matrix and the harmony of natural forms.

#### INTENDED USE

**Algeness®** is a reabsorbable filler indicated for the correction of tissue and/or dermal atrophies originating from trauma or pathological states. Treatment with **Algeness®** is part of reconstructive surgery, which more precisely aims to correct congenital malformations or functions compromised by trauma, disease, accidents, burns, wounds, often accompanied by loss of tissue substance. In addition, **Algeness®** is effective for the treatment of dermal atrophy and connective tissue deficits. **Algeness®** can also be used for deep filling of skin tissue. For **Algeness®** the clinical advantage is the reduction of the psychological impact, which can take on pathological consequences, caused by morphological changes in the patient/consumer's body.

Non-medical intended use: facial filling for the correction of skin imperfections.

In the aesthetic field, **Algeness®** is a resorbable filler indicated for the correction of skin imperfections caused by wrinkles, folds or sunken scars. In particular, the use of **Algeness VL** and **DF** is also specifically indicated for the correction of nasolabial wrinkles.

**Algeness® LD (low density - subdermal)** is indicated for repairs and superficial subdermal corrections surface repairs and subdermal corrections of hypotrophy or tissue lesions such as in the perioral area.

**Algeness® HD (medium density - subdermal)** is indicated for repairs and superficial subdermal corrections due to hypotrophy or tissue lesions in the perioral area.

**Algeness® VL (medium / high density - deep subdermal)** is indicated in the volumetric restoration for all forms of deep tissue hypotrophy with severe or modest loss of tone of the deep ligaments such as in the suborbital areas and in the anterior mandibular area, include correcting nasolabial wrinkles.

**Algeness® DF (high density - deep subdermal)** is indicated for all forms of volumetric restorations for all forms of severe deep tissue hypotrophy with severe loss of volume and tone of the deep ligaments such as in the suborbital areas and in the anterior mandibular area, include correcting nasolabial wrinkles.

#### INDICATIONS AND USE

**Algeness®** may only be used by doctors qualified in aesthetic medicine, dermatologists or doctors specialising in plastic, reconstructive and/or aesthetic surgery.

Before using the device, scan the QR code below to access training documentation on how to use the device, techniques and other useful information for safe use of **Algeness®**. The treatment must only be performed by competent doctors trained in the use of injection techniques.

Prior to use, connect the supplied backstop to one of the two syringes according to the relevant installation sheet Through fig.1 and fig.2.

Afterwards mix the gel back and forth between the two syringes at least ten times (10) for improving consistency and better results. Then, disconnect the empty syringe with the connector and attach the supplied needle to the pre-filled syringe.

**Algeness®** is a slowly absorbed filler and is totally biocompatible, it is indicated for the treatment of atrophic and hyperkeratotic changes of the skin and for the correction of wrinkles, skin folds, scars, or for problems associated with traumatic injuries.

Prior to the injection, an anamnestic examination should be performed to identify possible patient/consumer-related factors, such as current and/or previous treatments, procedures, accidents or conditions that may affect the procedure.

Before the injection, it is necessary to carry out a thorough cleansing and complete disinfection of the implant site. It is recommended to maintain a correct aseptis of the treated area before, during and after the insertion of the implant as well as to ensure a proper environment (clinic).

Place **Algeness®** under the site you want to correct at the level of the subdermis and at a variable depth depending on the agarose concentration. The amount to be injected into the subdermis and its depth is left to the judgement of the doctor. The amount of **Algeness**

applied at the same injection site during the session should not exceed 20 ml. **Algeness®** should not be used in excess (overcorrection and/or overdosage).

The correction can be performed with the usual techniques: linear threading technique, deep linear or fan technique. The injection must be carried out slowly for better positioning of the product in the desired areas with less trauma to the tissues.

The product is extruded by applying continuous and constant pressure on the plunger of the syringe while withdrawing the needle.

For an optimal correction, after implantation, always massage the treated area to make a homogeneous distribution of the injected material.

It is recommended to apply a cold compress to the treatment area to reduce any local undesirable reactions.

A monitoring period after administration is recommended to detect potential unwanted side effects.

For proper maintenance of the result, remind the patient/consumer of the importance of a follow-up a few months after treatment.

The patient/consumer must also be informed that the treatment can be suspended at any time at their request.

Inform the patient/consumer they must not apply cosmetics to the implant site for the next twelve hours and not to expose themselves to direct heat for the next few days (for example: exposure to the sun or to UVA and UVB rays, use of hair dryers or of hair dryers, heat radiated by fireplaces, saunas, etc.).

Do not exceed a maximum of 20 ml every 6 months.

For subsequent injections at the same site wait at least 15 days after the previous session - the same techniques must be followed.

If properly used, the product has a useful life of about 6 months and complete reabsorption in about 12 months.

#### CONTRAINDICATIONS

The use is contraindicated in all cases other than those specified in the INDICATIONS of the product, therefore, it must not be injected superficially at the level of the dermis or as a bolus.

Prior to the injection, assure that the patient/consumer does not experience particular hypersensitivity to any of the components in the product, otherwise not use on such subjects.

As with all fillers, the use of the product for the correction of particularly vascularized areas is not recommended, as it can increase the risk of compression and occlusion of the vessels and related phenomenology.

**Algeness®** must not be injected into blood vessels, as it could cause the onset of occlusion, local tissue necrosis or emboli.

**Algeness®** must not be injected in subjects with acute or chronic skin diseases in the areas to be corrected or in the immediate vicinity and, for precautionary reasons, in subjects with a history of anaphylactic reactions or severe allergies, in patients with severe organ or systemic diseases, including autoimmune disorders. The product is not intended for use persons under 18 years old, pregnant and or lactating women. It is recommended to avoid the combined use with other substances such as crosslinked fillers in the same treatment area.

The use of **Algeness®** by users not included in the category indicated in the previous paragraph, such as practitioners not trained, not qualified or accredited in the medical field or without a qualification in aesthetic medicine, dermatologists or doctors specialising outside the field of plastic, reconstructive and/or aesthetic surgery, is strictly prohibited.

#### PRECAUTIONS

As with all percutaneous procedures, the implant can entail risks of infection if the product is injected into anatomical areas where inflammatory or infectious processes are present, or without proper cleansing and disinfection of the area to be corrected.

In patient/consumers with bleeding and / or coagulation disorders or during treatments with anticoagulants, the product should be used with caution, as the injection may frequently cause local bleeding or bruising.

#### WARNINGS

Before the product is used, the patient/consumer must be provided with the instructions for use in order to be informed about the indications, professional use only, contraindications and possible side effects, including their treatment measures, that may occur with the implant.

Strictly follow the recommended injection techniques to minimise the likelihood of side effects.

Maintain proper hygiene during treatment to prevent infection.

If careful technique and massage is not performed following the injection, lumps may appear in the treated area.

Product is for SINGLE USE ONLY. Sterile product. Product must not be re-sterilized.

The stickers on the inner package, which indicate the sterility of the product should be red, otherwise DO NOT use the product.

In case of any irregularities with the packaging, isolate the non-compliant packaging and dispose of as special waste. Never use the product if it has not been stored correctly or if it has been frozen.

The unused material left in the syringe must be eliminated (disposed of) after treatment to avoid any risk of cross-infection due to use on other patient/consumers.

The doctor and / or the patient/consumer must report any serious incidents occurring in relation to the device to the manufacturer Ghimas and to the competent authority of the Member States where the user and / or patient/consumer is established.

FOR PROFESSIONAL USE ONLY.

DO NOT USE the product after the expiry date shown on the package.

Use the needle supplied or of the same specifications below.

Following the implantation, edematous reactions or skin redness, of a mild or moderate degree, rarely occur, however these will resolve completely within a few days.

Even if the application of **Algeness®** is not painful, the possibility of rare reports of transient pain due to the injection cannot be excluded, especially in particularly sensitive areas.

The use of resorbable and biocompatible fillers may entail some side effects. Common side effects and recommended treatments are listed below.

- **Nodules:** Nodules may appear within the first four weeks after injection and are single, well-confined and non-inflammatory. Prolonged massage or dispersal of the nodule through a cannula accelerates its biodegradation.
- **Hardening:** the slow and precise injection technique can prevent this effect. If it occurs, it is necessary to perform a gentle massage, and if necessary apply warm compresses. In some cases it is advisable to proceed with a broad-spectrum antibiotic treatment, for which medical advice should be sought.
- **Inflammation:** Apply cold compresses and use anti-inflammatory creams if necessary. The inflammation should subside within a few days.
- **Immune response:** Agarose gel is a pure polysaccharide and does not contain any proteins, cross-linking agents or chemicals; therefore, the immune reaction is extremely low but can be managed by a doctor's antihistamine treatment if necessary.
- **Swelling:** any trauma, including injections, can cause temporary swelling. Applying cold and sleeping with the head slightly elevated can promote fluid drainage and reduce swelling. In some cases, gentle massage or the use of anti-inflammatory creams/gel is recommended.
- **Pain:** Typically disappears immediately after the injection act, but in other cases can be relieved with painkillers such as paracetamol or ibuprofen. Avoid aspirin and other anticoagulants to reduce the risk of haematoma.
- **Haematoma:** Bruising or haematoma formation is common. Apply cold compresses and use arnica-based creams to accelerate bruise resolution. Avoid sun exposure to affected areas.
- **Edema:** Apply cold compresses and drink plenty of water to facilitate lymphatic drainage. Reduce salt intake to avoid further water retention.
- **Overdosage:** As **Algeness®** is a hydrocolloid, there is no possibility of an unexpected increase in volume. However, if it is certain that there is an overcorrection after injection, the filler can be dispersed with an injection of warm saline and massage in the initial period.

The doctor must be consulted immediately if necessary and if adverse reactions and side effects occur.

In addition to the common side effects, other less common effects may also occur, which include but are not limited to:

- Capsule formation or capsular contractures,
- Infection,
- Superficial wound,
- Hyperpigmentation,
- Seroma,
- Increased pressure within the compartment,
- Granuloma.

These side effects are less likely but should be monitored closely. If any of these symptoms or others occur, it is essential to consult a doctor immediately.

#### PACKAGING

The product is supplied in a pouch with a special double syringe Luer Lock "total emptying". Each box contains two pouches and four sterile needles.

The volume contained in the syringe is indicated on the pouch label.

The volume on the syringe is purely indicative.

#### STORAGE

Conserve the product at room temperature, between 6°C and 30°C. If necessary, the product may be stored in a refrigerator, the product must not be frozen.

#### HARMONISED STANDARDS AND COMMON SPECIFICATIONS

**Algeness®** has been developed and manufactured in accordance with Regulation (EU) 2017/745 as amended and according to the Common Specification EU Regulation 2022/2346 as amended.

**Algeness®** is also compliant with the following applicable harmonised standards updated to 6 March 2024 (EU Implementing Decision 2024/815):

<b>EN ISO 13485</b> "Medical devices - Quality management systems - Requirements for regulatory purposes"
<b>EN ISO 14971</b> "Medical devices - Application of risk management to medical devices"
<b>EN ISO 15223-1</b> "Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements"
<b>EN ISO 10993-10</b> "Biological evaluation of medical devices Part 10: Tests for irritation and skin sensitization"
<b>EN ISO 10993-12</b> "Biological evaluation of medical devices Part 12: Sample preparation and reference materials"
<b>EN ISO 10993-17</b> "Biological evaluation of medical devices Part 17: Establishment of allowable limits for leachable substances"
<b>EN ISO 10993-18</b> "Biological evaluation of medical devices Part 18: Chemical characterization of materials"
<b>EN ISO 10993-23</b> "Biological evaluation of medical devices Part 23: Tests for irritation"
<b>EN ISO 11137-1</b> "Sterilization of health care products — Radiation — Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices"
<b>EN ISO 11137-2</b> "Sterilization of health care products — Radiation — Part 2: Establishing the sterilization dose"
<b>EN ISO 11737-1</b> "Sterilization of health care products — Microbiological methods — Part 1: Determination of a population of microorganisms on products"
<b>EN ISO 11737-2</b> "Sterilization of health care products — Microbiological methods — Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process"
<b>EN ISO 11607-1</b> "Packaging for terminally sterilized medical devices — Part 1: Requirements for materials, sterile barrier systems and packaging systems"
<b>EN ISO 11607-2</b> "Packaging for terminally sterilized medical devices — Part 2: Validation requirements for forming, sealing and assembly processes"

On the basis of MDR 2017/745 and subsequent modifications, the physician must ensure the traceability of the product, i.e. he must be able to trace the patient/consumer to whom he applied the product. It is recommended to attach the product label to the patient/consumer's medical record.

Consult the attachment of the instructions for completing the implant card, take the inserted card and deliver it **completed** to the patient/consumer.

Fill in the specific section indicated in the consumer instructions and return it to the consumer.

It is recommended to carefully read the Summary relating to safety and clinical performance of **Algeness®** published on the Ghimas website, [www.ghimas.com](http://www.ghimas.com), in the dedicated area and/or on the European database "Eudamed".

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Ghimas S.p.A.
Advanced Aesthetic Technologies Inc. – Brookline MA, USA